REAL-TIME IMMUNOASSAYS

Rapid I Precise I Simple



Passage du Cardinal 13b 1700 Fribourg Switzerland

info@elionova.com www.elionova.com



ト ス 田 ス の の の

- About us
- Our Products
- Our Service
- Our Technology
- Evaluation approaches
- Our Workflow
- Assay performance
- Assay matrices
- Our Applications
- Testimonials



Precise & Rapid diagnostics for Research Use, IVD & Beyond

About Us

ElioDX™ wash-free technology delivers real-time, laboratory-quality, immunoassays in just 10 minutes through a simple, one-step, disposable system.

Designed for unmatched precision, speed, and versatility, ElioDXTM supports a wide array of applications across Research & Development (R&D), In-vitro Diagnostics (IVD), Food Safety, Agriculture, Drug Development and Quality Control.



Why Choose Us?



SPEED

10 min to result

10x faster than the industry benchmark



HIGH PRECISION and ACCURACY

Quantitative immunoassays.

Real-time results provides further information on reaction.



SIMPLICITY

Single pipetting step, **no washing or incubation.**

Cuvettes contain all reagents in dried form



FLEXIBILITY

Offers users the possibility to develop own assays.

Bring **existing** assays to ElioDX platform.



DIGITAL RECORDS

Digital platform for assay setup, evaluation and data storage.

Results
accessible from
any computer
and smart phone.

OUR PRODUCTS


ElioDXTM INSTRUMENT

Compact & Portable
Laser Evanescences technology
Photon counting
Highly precise mechanics



ElioDXTM DISPOSABLES

ElioDXTM PRE-MADE ASSAYS



the same time



Human IgG I-FABP FSH IL-6 PCT IL-27 CRP CD-64 SAA MPO

+ 50 applications a year





single step





Our DIY kit,
NeutrAvidin kit and
Protein G kit provide a
fully customisable
workflow—letting users
choose their own
workflow to suit their
specific application.







OUR SERVICE

Custom Assays

Elionova offers custom assay development tailored to your research needs. From initial request to validated delivery, our team ensures each assay is optimized for accuracy, speed, and seamless integration with your workflow.



Client Request Initiation

Expert feasibility review


Proposal and Assay setup review

Assay development and Validation

Assay delivery and Support

Conjugation


Elionova offers a reliable protein conjugation service to Allophycocyanin (APC), Cy5, or Alexa Fluor 647. From planning to purification and quality control, we ensure high-performance conjugates tailored to your application.





A clever combination of:

Immunoassay methods
Laser Evanescence Technology
Photon Counting
Highly Precise Mechanics



ELIONOVA

KEY POINTS:


- **Evanescent laser field** excites molecules only at the bottom of the well, where the immune reaction takes place
- Specific signal increase over time exceeds the background signal without analyte
- No washing steps required

EVALUATION APPROACHES

Dynamics (1) vs steady state (2)

ElioDX™ is a versatile and flexible system that enables:

- (1) dynamic monitoring fluoresent signal overtime (10 min) and plot slope vs concentration for a kinetic standard curve;
- (2) steady state capturing endpoint signal in just 2 minutes and plot fluorescent count vs. concentration for traditional standard curve perfect for high-throughput or simplified protocols.




The measured **slope** (dynamic) or **count** (steady state) of the sample of interest is referenced against the standard curve to calculate the analyte concentration.

OUR WORKFLOW

How to operate a pre-made **ELIODX**TM **IMMUNOASSAY**


Each ElioDXTM cuvette includes a **unique QR code** for easy registration and metadata entry, such as sample IDs and sample matrix details. Samples are pipetted directly into pre-loaded cuvettes—no incubation or washing required. The ElioDXTM reader then performs **10 min real-time**, quantitative fluorescence measurements. Results are automatically uploaded to **the ElioDXTM cloud**, where they're processed and available in Excel or CSV format.



ASSAY PERFORMANCE


The ElioDX[™] system demonstrates analytical sensitivity comparable to, or exceeding other immunoassays, such as ELISA and chemiluminescence, with high efficiency from its no-wash protocol.

ElioDXTM Human Follicle Stimulating Hormone (FSH)




Development of fluorescent signal over time

(0-10 min) for different concentrations of FSH




Analytical limit of detection is 0.3 U/L (clinical FSH cut-off is 1 U/L)




correlate well
(R² = 0.99) to the values
obtained using **Roche COBAS** analyser.

ElioDXTM Human Procalcitonin (PCT)




Development of fluorescent signal over time

(0-10 min) for different concentrations of PCT



Analytical limit of detection is 0.1 ng/mL (clinical PCT cut-off is 0.1 ng/ml)




Elionova PCT results

correlate well (R² = 0.95) to the values obtained using ECLIA analyser.

ASSAY MATRICES

The ElioDX™ system is compatible with a wide range of biological sample types, including blood, saliva, urine, cerebrospinal fluid, and cell culture supernatants.

The **example** below shows different sample matrices—**buffer, heparin plasma,** and **serum**—tested in our **ElioDXTM Human PCT assay**, with the signal development plot on the left and the corresponding standard curve on the right for each matrix.



OUR APPLICATIONS

The ElioDX™ system detects a broad range of analytes—including antibodies, proteins, nucleotides, and drugs—using flexible assay formats like direct, indirect, sandwich, and competitive, similar to ELISA.

Analytes







Protein



Nucleotide



Assay formats



Direct



Indirect



Sandwich



Competitive

Protein interaction characterization

Our technology provides functionality comparable to surface plasmon resonance (SPR) analysis, but with exceptionally simple handling and low cost.





Kinetics How fast does binding occur?



Concentration How much of the molecule is present?



Specificity How selective is the interaction?

Application fields



RESEARCH

Cell biology Infections **Immunology** Immuno-oncology Neurobiology



DIAGNOSTICS

Veterinary IVD Human IVD Food safety Agriculture



PHARM

Drug discovery Development **Quality control**

•••



TESTIMONIALS



Dr. Anthony Croxatto,
Director of the Department
of ADMED Medical Microbiology,
Swiss Reference Centre
for Tick-borne diseases.

"Elionova's real time immunoassay technology has dramatically improved our immunoassays for detection of tick-borne diseases."



Peter Spies,
Research Scientist
University of Applied Sciences and
Arts Northwestern Switzerland

"We anticipate a massive gain in time to result, material as well as hands-on time making the ElioDX™ technology much more economic than our traditional methods like ELISA."

