

BioConcept Standard Agarose Type LE for Routine Gel Electrophoresis

Cat. No.	Product	Size	Offer Price CHF
7-01P02-R	BioConcept Standard Agarose Type LE for routine gel electrophoresis	500g	270.00

This promotion is valid until December 31, 2024. No other discounts will apply. All prices in CHF excl. VAT. Please mention quotation number **BAP1224** when ordering. If you are interested in another package size of our agarose (25g, 100g, or 1000g) please ask us for a quotation.

TBE, TAE and SSC buffers and Water for Molecular Biology

Cat. No.	Product	Size	List Price CHF
3-07F01-I	10x TBE buffer	500ml	16.80
3-07F03-I	50x TAE buffer	500ml	24.00
3-07F00-I	20x SSC buffer, pH 7.0	500ml	14.70
3-07F04-H	Water for Molecular Biology	100ml	9.00
3-07F04-I	Water for Molecular Biology	500ml	16.00


All prices in CHF excl. VAT.



NEB's DNA Markers and Ladders – built for perfect migration!
Please check reverse side for available options and a free sample.



NEB's Markers and Ladders – built for perfect migration!




Ready-to-use, pre-diluted formats for added convenience:

- **Quick-Load® Purple:** containing purple loading dye
- **Quick-Load®:** containing Bromophenol Blue
- **TriDye™:** containing three dyes for better visualization
- **For Safe Stains:** compatible with GelRed®, GelGreen® and SYBR® precast gels

Included with NEB DNA Markers & Ladders: Gel Loading Dye, Purple (6X)

- Brighter, sharper bands vs. conventional loading buffers
- No UV shadow allows for publication-grade images
- Stops enzyme reactions



Gel Loading Dye, Purple (6X), no SDS, is included in most DNA ladders and markers. This unique loading buffer does not cast a UV shadow over the underlying bands, unlike conventional Bromophenol Blue containing loading buffers.

Ordering Information

PRODUCT	NEB #	SIZE
1 kb Plus:		
1 kb Plus DNA Ladder	N3200S/L	200/1,000 gel lanes
Quick-Load Purple 1 kb Plus DNA Ladder	N0550S/L	250/750 gel lanes
Quick-Load 1 kb Plus DNA Ladder	N0469S	250 gel lanes
TriDye 1 kb Plus DNA Ladder	N3270S	250 gel lanes
1 kb Plus DNA Ladder for Safe Stains	N0559S	250 gel lanes
1 kb:		
1 kb DNA Ladder	N3232S/L	200/1,000 gel lanes
Quick-Load Purple 1 kb DNA Ladder	N0552S/L	125/375 gel lanes
Quick-Load 1 kb DNA Ladder	N0468S/L	125/375 gel lanes
TriDye 1 kb DNA Ladder	N3272S	125 gel lanes
Quick-Load 1 kb Extend DNA Ladder		
Quick-Load 1 kb Extend DNA Ladder	N3239S	125 gel lanes
100 bp:		
100 bp DNA Ladder	N3231S/L	100/500 gel lanes
Quick-Load Purple 100 bp DNA Ladder	N0551S/L	125/375 gel lanes
Quick-Load 100 bp DNA Ladder	N0467S/L	125/375 gel lanes
TriDye 100 bp DNA Ladder	N3271S	125 gel lanes
50 bp:		
50 bp DNA Ladder	N3236S/L	200/1,000 gel lanes
Quick-Load Purple 50 bp DNA Ladder	N0556S	250 gel lanes
Low Molecular Weight:		
Low Molecular Weight DNA Ladder	N3233S/L	100/500 gel lanes
Quick-Load Purple Low Molecular Weight DNA Ladder	N0557S	125 gel lanes
Fast DNA Ladder:		
Fast DNA Ladder	N3238S	50 gel lanes



Order your free sample of Quick-Load Purple 1 kb Plus DNA Ladder (10-20 gel lanes) #N0550G at BioConcept!