





Thermal Cycler Compatibility Guide



	96-Well Semi-Skirt 0.2mL	96-Well Full-Skirt 0.1mL		384-Well Full-Skirt A24 Notch		384-Well Full-Skirt A24, P24 Notches		
	Clear Wells	Clear Wells	White Wells	Clear Wells	White Wells	Clear Wells	White Wells	
APPLIED BIOSYSTEMS® (ABI)								
3100, 3130XL, 3500, 3500XL, 3730, 3730XL				●	●			
GeneAmp 2700, 2720	●							
GeneAmp 9600, 9700	●							
MiniAmp, MiniAmp Plus	●							
ProFlex	●							
Proflex, Multiblock system				●	●			
QuantStudio 3, 5, 6, 7, 12 K Real-time PCR System (with 0.2mL blocks)	●			●	●			
SimpliAmp	●							
Veriti 384-well				●	●			
Veriti 96-well (with 0.2mL block)	●							
ViiA7, 7900 HT FAST				●	●			
AGILENT®								
AriaMx qPCR System		●	●					
Mx3000P, Mx3005P, Mx4000	●							
Robocycler Gradient	●	●	●					
Surecycler 8800	●			●	●			
ANALYTIK JENA®								
Trobot 96, SpeedCycler2 SP, SPR								
Alpha SC								
Flexcycler2	●	●	●	●	●			
T1 Thermocycler, T Gradient, T One, T Advanced, T Professional series	●	●	●	●	●			
qTOWER3 (G, touch), Toptical	●			●	●			
BIOER®/ BULLDOG®								
GeneExplorer	●							
GeneTouch 384				●	●			
GeneTouch, GeneTouch Plus	●							
LifeECO	●							
LifeTouch	●							
LifePro	●							
BIOMETRA®								
Biometra T Advanced				●	●			
Biometra T Advanced 96, 96 G, 96 S, 96 SG	●	●	●					
Biometra T One	●	●	●					
Biometra T Professional series	●	●	●					
Biometra Uno	●	●	●					




Thermal Cycler Compatibility Guide



	96-Well Semi-Skirt 0.2mL	96-Well Full-Skirt 0.1mL		384-Well Full-Skirt A24 Notch		384-Well Full-Skirt A24, P24 Notches		
	Clear Wells	Clear Wells	White Wells	Clear Wells	White Wells	Clear Wells	White Wells	
BIO-RAD® / MJ RESEARCH®								
C1000, C1000 Touch, S1000	●	●	●	●	●			
CFX Connect, CFX96, CFX96 Touch		●	●					
CFX384, CFX384 Touch				●	●			
CFX96 Touch Deep Well	●	●	●					
Chromo4	●	●	●					
DNA Engine Opticon and Opticon 2		●	●					
DNA Engine, DNA Engine Tetrad, Tetrad 2		●	●	●	●			
DNA Engine Dyad, Dyad Disciple		●	●	●	●			
iCycler	●							
iQ, iQ-4, iQ-5, myiQ PCR System	●							
MJ Research PTC200, 220, 221, 225, 240	●			●	●			
MyCycler	●							
PTC 200		●	●	●	●			
T100	●							
CORBETT® / QIAGEN®								
Corbett Research (Qiagen) Palm Cycler	●							
Corbett Research (Qiagen) Palm Cycler 384				●	●			
EPPENDORF®								
MasterCycler ep 384 / Pro 384				●	●			
Mastercycler ep realplex	●	●	●					
MasterCycler ep, ep gradient, Pro, Pro S	●	●	●					
Nexus, Nexus gradient, Nexus SX1, Nexus GSX1	●	●	●					
GE / MEGABACE®								
GE Healthcare Amersham MegaBACE 4000				●	●			
GE Healthcare Amersham MegaBACE 500, MegaBACE 1000 mark 2		●	●					
ROCHE®								
LightCycler 480 (384 well)						●	●	
TECHNE® / COLE-PARMER®								
Flexigene	●	●	●	●	●			
Techne Genius, Genius Quad	●	●	●	●	●			
Techne PCRmax Alpha cycler 1, 2, 4				●	●			
Techne Prime, PrimeG, Prime Elite, Prime Elite Satellite	●			●	●			
Techne Quantica		●	●	●	●			
Techne TC412, TC512		●	●	●	●			
Techne TC-PLUS	●	●	●					
Touchgene, Touchgene Gradient	●	●	●					



Thermal Cycler Compatibility Guide



	96-Well Semi-Skirt 0.2mL	96-Well Full-Skirt 0.1mL		384-Well Full-Skirt A24 Notch		384-Well Full-Skirt A24, P24 Notches	
	Clear Wells	Clear Wells	White Wells	Clear Wells	White Wells	Clear Wells	White Wells
OTHER BRANDS / MODELS							
Apollo ATC-201	●						
BIONEER AllInOneCycler	●						
DLab TC1000-G	●	●	●				
G-Storm GS1, GS2, GS4	●	●	●				
HiMedia® Prima-96	●						
Labnet MultiGene, OptiMax							
Idaho Technology LightScanner		●	●	●	●		
peqLab peqSTAR 384X							
peqLab peqSTAR 96X	●						
Sensoquest Labcycler	●			●	●		
Takara Dice touch, Gradient	●						
ThermoFisher Hybaid MBS, Px2/PxE	●						
VWR Collection UNO96	●	●	●				
VWR peqlab peqSTAR 96X		●	●				