

Fall Winter Promotion 2026

High-Fidelity PCR with NEB

Enjoy now an exclusive 25% discount on your
purchase of Q5 Polymerases!



Provided by BioConcept AG

Q5 DNA Polymerase

Benefit from the quality of Q5 DNA Polymerase

Q5 High-Fidelity DNA Polymerase sets the standard for both fidelity and robust performance. With fidelity >280 times higher than Taq DNA Polymerase, Q5 results in ultra-low error rates. Q5 is a designed polymerase, improving speed, fidelity, and reliability of performance.



Extremely low error rates

At ~280X higher than Taq, Q5 offers unparalleled fidelity for your most important samples.



Robust amplification with minimal optimization

High specificity and yield are absolute requirements for today's molecular biology techniques.



Superior coverage, regardless of GC content

Q5 displays superior performance for a wide range of high-GC or high-AT templates.



Shorter PCR protocols

Incorporating the SSo7d processivity-enhancing domain enables shorter extension times, as low as 5 seconds per kb.




Room Temperature setup

Aptamer-based hot start requires no initial denaturation step and enables room temperature setup.



Universal annealing temperature

Compatible with a universal annealing temperature protocol.



Amplification of long templates

Reliably amplify simple templates up to 30 kb and complex templates up to 20 kb.



Giving Q5 a try?

Request a **free sample** from BioConcept:
info@bioconcept.ch

New England Biolabs — Fall Winter Promotion 2026

Available Q5 DNA Polymerases

Cat. No	Q5 DNA Polymerase	List Price	Offer
M0491S	Q5 High-Fidelity DNA Polymerase, 100 units	161.00	120.75
M0491L	Q5 High-Fidelity DNA Polymerase, 500 units	707.00	530.25
M0493S	Q5 Hot Start High-Fidelity DNA Polymerase, 100 units	190.00	142.50
M0493L	Q5 Hot Start High-Fidelity DNA Polymerase, 500 units	864.00	648.00
M0492S	Q5 High-Fidelity 2X Master Mix, 100 rxns (50 µl vol)	265.00	198.75
M0492L	Q5 High-Fidelity 2X Master Mix, 500 rxns (50 µl vol)	1060.00	795.00
M0494S	Q5 Hot Start High-Fidelity 2x Master Mix 100 rxns (2x1.25ml)	298.00	223.50
M0494L	Q5 Hot Start High-Fidelity 2x Master Mix, 500 rxns (10x1.25ml)	1192.00	894.00
M0494X	Q5 Hot Start High-Fidelity 2x Master Mix, 500 rxns (1x12.5ml)	1192.00	894.00
E0555L	Q5 High-Fidelity PCR Kit, 200 rxns	530.00	397.50
E0555S	Q5 High-Fidelity PCR Kit, 50 rxns	156.00	117.00
M0515S	Q5U Hot Start High-Fidelity DNA Polymerase, 100 units	255.00	191.25
M0515L	Q5U Hot Start High-Fidelity DNA Polymerase, 500 units	1081.00	810.75
M0500S	Q5 Blood Direct 2x Master Mix, 100 rxns (1x1ml)	170.00	127.50
M0500L	Q5 Blood Direct 2x Master Mix, 500 rxns (5x1ml)	740.00	555.00
M0578L	Q5 Quick-Load High-Fidelity 2x Master Mix, 500 rxns (10x1.25ml)	1060.00	795.00
M0580L	Q5 Quick-Load Hot Start High-Fidelity 2x Master Mix, 500 rxns (10x1.25ml)	1192.00	894.00
M2499S	Q5-XT Hot Start High-Fidelity 2x Master Mix, 100 rxns (2x1.25ml)	308.00	231.00
M2499L	Q5-XT Hot Start High-Fidelity 2x Master Mix, 500 rxns (10x1.25ml)	1220.00	915.00
M0544S	NEBNext Ultra II Q5 Master Mix, 50 rxns	155.00	116.25
M0544L	NEBNext Ultra II Q5 Master Mix, 250 rxns	615.00	461.25
M0544X	NEBNext Ultra II Q5 Master Mix, 500 rxns	1169.00	876.75
M0597S	NEBNext Q5U Master Mix, 50 rxns	193.00	144.75
M0597L	NEBNext Q5U Master Mix, 250 rxns	766.00	574.50

Prices in CHF

**Save 25% on
Q5 Polymerases**

Simply enter the promo code "Q52026"
when ordering at BioConcept.

*The 25% discount is applicable to the list price of Q5 products listed above. It is valid until December 31, 2026, and cannot be combined with other discounts. For more details, please contact [BioConcept](mailto:info@bioconcept.ch).

BioConcept Ltd.
Paradiesrain 14
Postfach 427
4123 Allschwil
Switzerland

Tel. +41 (0)61 486 80 80
info@bioconcept.ch
www.bioconcept.ch

