

New Products from CST®: September 2025


We are continually expanding our portfolio to meet your research needs. Check out the new products released this month!

InTraSeq™ 3' Conjugate Antibody Cocktail 2 #83302

The second InTraSeq™ 3' Conjugate Antibody Cocktail!

A new antibody cocktail for InTraSeq (Intracellular Protein and Transcriptomic Sequencing), which enables signal pathway analysis at the single-cell level, has been released! Developed by CST, InTraSeq single-cell analysis enables identification of signaling pathways involved in disease development and elucidation of molecular mechanisms in just one experiment. By simultaneously detecting both intracellular and cell surface proteins as well as RNA levels in thousands of cells, signal pathway analysis using the transcriptome at the single-cell level is possible.

KEYWORDS: single-cell analysis, signal pathway analysis, post-translational modification




SCA-ITS3: A-204 cells, 293 cells, serum-starved (ss) HT-1080 cells, and ss HT-1080 cells treated with Human Transforming Growth Factor β3 (hTGF-β3) #8425 (100 ng/mL, 30 min) were processed in accordance with the InTraSeq™ 3' protocol and analyzed using the Seurat computational package. The leftmost figures in each row show the Uniform Manifold Approximation and Projection (UMAP) of the annotated cells. The top right FeaturePlot displays the Phospho-SMAD2 (Ser465/Ser467) (E8F3R) Rabbit mAb (InTraSeq™ 3' Conjugate 3061) #34360 expression. The bottom figures show a subclustering analysis where ss HT-1080 cells and ss HT-1080 + hTGF-β3 cells were computationally selected to show expression of Phospho-SMAD2 (Ser465/Ser467) (E8F3R) Rabbit mAb (InTraSeq™ 3' Conjugate 3061) #34360 with ss HT-1080 cells in the middle figure and ss HT-1080 + hTGF-β3 cells in the right figure. To generate similar plots, use the FeaturePlot command template: `FeaturePlot(your_object, features = your_protein, split.by = "sample", max.cutoff = "q95", min.cutoff = "q5", cols = c("grey90", "#b22222"))`.

CXCR4 (F1K2M) Rabbit mAb #66825

Chemokine receptor involved in promoting cancer cell proliferation, survival, and motility

CXCR4 is a chemokine receptor belonging to the G protein-coupled receptor family. It is activated by CXCL12, also known as SDF1. CXCR4 mediates the homing of progenitor cells in the bone marrow and their recruitment to sites of injury, as well as regulating the migration, proliferation, and survival of cancer cells. It is also being investigated as a potential therapeutic target in autoimmune diseases and cancer. #66825 is suitable for WB and IHC-P using mouse samples.

KEYWORDS: chemokine, cancer microenvironment, metastasis



IHC-P: IHC analysis of paraffin-embedded human gastric adenocarcinoma using #68008. IHC analysis of paraffin-embedded A20 syngeneic tumor using #66825.



Learn more at
cellsignal.com/newproducts

PRODUCT		APPLICATIONS	REACTIVITY
Adhesion/Extracellular Matrix			
34720S	ADAM10 (F9A4J) Rabbit mAb	IHC-P	H
34720T			
67360S	COL17A1 (F5U1Z) Rabbit mAb	WB, IHC-P	H
67360T			
19775SF	DPP4/CD26 (D6D8K) Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H
77093S	GPNMB (E7U1Z) Feline Chimeric mAb	IF-F	M
77093T			
49553S	GPNMB (E7U1Z) Mouse Chimeric mAb	IF-F	M
49553T			
23275S	Mesothelin (D9R5G) XP® Rabbit mAb (Alexa Fluor® 488 Conjugate)	IHC-P	H
50141S	Mesothelin (D9R5G) XP® Rabbit mAb (Alexa Fluor® 647 Conjugate)	IHC-P	H
91115SF	Phospho-Paxillin (Tyr118) (E9U9F) Rabbit mAb (BSA and Azide Free)	WB, IF-F, IF-IC	H, M, R, Mk, (C)
96408SF	Periostin (F7U5N) Rabbit mAb (BSA and Azide Free)	WB, IHC-P	M
86662S	Thy1/CD90 (D3V8A) Rabbit mAb (Alexa Fluor® 647 Conjugate)	IHC-P	H, M, R
52393S	VE-Cadherin (E6N7A) Rabbit mAb (Alexa Fluor® 488 Conjugate)	IHC-P	H
19948S	VE-Cadherin (E6N7A) Rabbit mAb (Alexa Fluor® 647 Conjugate)	IHC-P	H
Alzheimer's disease			
66973S	CD11b/ITGAM (E6E1M) Feline Chimeric mAb	IF-F	M
66973T			
78923S	CD11b/ITGAM (E6E1M) Horse Chimeric mAb	IF-F	M
78923T			
90153S	CD11b/ITGAM (E6E1M) Mouse Chimeric mAb	IF-F	M
90153T			
81911S	P2Y12 (E9J1J) Mouse Chimeric mAb	IF-F	M
81911T			
31944S	TREM2 (E7O7Q) Rabbit mAb (Alexa Fluor® 488 Conjugate)	FC-L	M
44724S	TREM2 (E7O7Q) Rabbit mAb (Alexa Fluor® 647 Conjugate)	FC-L	M
Apoptosis			
50433SF	Phospho-RIP3 (Thr231/Ser232) (E7S1R) Rabbit mAb (BSA and Azide Free)	WB, IF-IC	M
Autoimmunity and Inflammation			
59098S	IL-1RA (F4K3K) Rabbit mAb	WB	H
59098T			
17382SF	p47phox (D21F6) Rabbit mAb (BSA and Azide Free)	WB, IHC-P	H, (Mk)
Autophagy			
78080SF	Phospho-mTOR (Ser2448) (D9C2) XP® Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H, M, R, Mk, (C, Pg, Hr)
44380SF	Mucopolin-1 (F8F9Q) Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H
24647SF	Phospho-SQSTM1/p62 (Ser349) (E7M1A) Rabbit mAb (BSA and Azide Free)	WB, IHC-P, IF-IC, FC-FP	H, M
89254T	TFEB/TFE3 Signaling Antibody Sampler Kit	-	-
Biochemical Metabolism and Signaling			
35190S	AFP (F2E9D) Rabbit mAb	WB, IHC-P	H
35190T			
43596S	AFP (F8P8E) Rabbit mAb	IHC-P	H
43596T			
63435SF	ATF-4 (E4Q4E) Mouse mAb (BSA and Azide Free)	WB, IF-IC, FC-FP	H
89633SF	ATGL (30A4) Rabbit mAb (BSA and Azide Free)	WB, IHC-P, IF-IC	M
16991S	ATP1B3 (F1R8N) Rabbit mAb	WB, IHC-P	H
16991T			
39539SF	Phospho-CAD (Ser1859) (D5K5W) Rabbit mAb (BSA and Azide Free)	WB, IHC-P, SEC	H, M, R
78571S	FATP1/SLC27A1 (F1Y1E) Rabbit mAb	WB	H
78571T			
62511T	Fatty Acid Binding Protein Antibody Sampler Kit	-	-
82536S	GSTT1 (F6G2R) Rabbit mAb	WB	H, M, R
82536T			
86051S	Iduronate 2-Sulfatase (F5F9E) Rabbit mAb	WB	H, M, R
86051T			
86188SF	INPP4b (D9K1B) XP® Rabbit mAb (BSA and Azide Free)	WB, IHC-P	H
94645S	MLH1 (F1W3E) Rabbit mAb	WB	H, M
94645T			
85653T	Pancreatic Development Antibody Sampler Kit	-	-
96703S	SGPP1 (F2B6W) Rabbit mAb	WB, IP, FC-FP	H
96703T			
32368S	TYRP1 (F7P2J) Rabbit mAb	WB, IHC-P, IF-F, IF-IC	H, M
32368T			
Ca, cAMP & Lipid Signaling			
85303S	IP3R3/ITPR3 (F3L2A) Rabbit mAb	WB	H, M, R
85303T			
32450S	Phospho-PKCζ (Thr410)/λ (Thr412) (F2E5C) Rabbit mAb	WB, IP	H, M, R, Mk
32450T			
16685SF	Phospho-PLCγ1 (Ser1248) (D25A9) Rabbit mAb (BSA and Azide Free)	WB, IHC-P, IF-IC, FC-FP	H, M, Mk, (R)



Learn more at
cellsignal.com/newproducts



Cell Signaling
 TECHNOLOGY®

PRODUCT	APPLICATIONS	REACTIVITY
Cell Autonomous Innate Immunity		
62664SF IRF-3 (D6I4C) XP® Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H, Mk
37230SF STING (E9X7F) Rabbit mAb (BSA and Azide Free)	IF-F, IF-IC, FC-FP	H, M
Cell Cycle		
93597S CDK4 (F8L5F) Rabbit mAb	WB, IHC-P	H, M, R
93597T		
62382SF CDK6 (D4S8S) Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H, Mk
33667C PathScan® RP CDK6 Sandwich ELISA Kit	ELISA	H, Mk
25211SF CDT1 (D10F11) Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H, Mk
79338S Geminin (F9T4M) Rabbit mAb	WB, IP	M
79338T		
75886SF p16 INK4A (F9T1L) Rabbit mAb (BSA and Azide Free)	WB, IHC-P, IF-IC, FC-FP	H
34418S PrimPol (F5E9I) Rabbit mAb	WB, IP	H
34418T		
69017SF Phospho-RPA32/RPA2 (Ser8) (E5A2F) Rabbit mAb (BSA and Azide Free)	WB, IF-IC, FC-FP	H
89040S TERT (F3C6F) Rabbit mAb	WB	H
89040T		
Chromatin Acylation and Phosphorylation		
89770SF SRC-3 (D1F11) Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H, M, Mk
Chromatin Methylation and Ubiquitylation		
51519SF APOBEC3A (UMN-13) Rabbit mAb (BSA and Azide Free)	WB, IHC-P	H
Cytoskeletal		
40399SF MYH11 (E4Q4L) Rabbit mAb (BSA and Azide Free)	WB, IF-F, IF-IC	H, M, R
95377SF Phospho-Myosin Light Chain 2 (Ser19) Mouse mAb (BSA and Azide Free)	WB, IF-IC	H, M, R, B, Pg
18751S Phospho-NCK1 (Ser85) (F8PIY) Rabbit mAb	WB, IP, FC-FP	H, M
18751T		
DIFFERENTIATION - Transcription Factor		
27425SF GATA-1 (D24E4) XP® Rabbit mAb (BSA and Azide Free)	WB, IF-IC, FC-FP	H
57731S KLF2 (E7K8Y) Rabbit mAb (Alexa Fluor® 647 Conjugate)	FC-FP	H
46846SF Sox2 (D9B8N) Rabbit mAb (BSA and Azide Free)	WB, IF-F, IF-IC, FC-FP	H, M
DNA Damage/Senescence		
51030SF Phospho-Rad50 (Ser635) (F5H8B) Rabbit mAb (BSA and Azide Free)	WB, FC-FP	H, R
54069S Rad51C (F6D1F) Rabbit mAb	WB, IP	H
54069T		
77441SF DeltaN p63 (E6Q3O) Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H
51201SF WRN (8H3) Mouse mAb (BSA and Azide Free)	WB, IF-IC	H, M
General Transcription Factors		
48494SF HEXIM1 (D5Y5K) Rabbit mAb (BSA and Azide Free)	WB, IF-IC, FC-FP	H, Mk
96999S INTS11 (F2P6Z) Rabbit mAb	WB, IP	H, M, R, Mk
96999T		
76072S NPAT (F1D9Q) Rabbit mAb	WB, IP, IF-IC	H, Mk
76072T		
Growth factors & Cytokines		
57408SF HGF β (D6S7D) XP® Rabbit mAb (BSA and Azide Free)	WB, IHC-P	H
HIPPO		
24663S MST2 (F4C3I) Rabbit mAb	WB, IP	H, M, R, Mk
24663T		
Immune Checkpoints and Ligands		
26043S B7-H3 (F9C1J) Rabbit mAb	WB	H, M, R
26043T		
77017S CD27 (F3R2V) Rabbit mAb (Alexa Fluor® 488 Conjugate)	IF-F, IF-IC, FC-L	M
56043S		
19040S CD40 (E2Z7J) Rabbit mAb (Alexa Fluor® 488 Conjugate)	IHC-P	M
Inflammasome		
81980S IFI16 (F6T5U) Rabbit mAb	WB, IP, IF-IC	H
81980T		
58464S NLRP7 (F9U8L) Rabbit mAb	WB, IP	H
58464T		
Lymphocyte Signaling		
71051P BLK (E8T1B) Rabbit mAb (Biotinylated)	E-P	H
79809S CD79A (F3Z1H) Rabbit mAb	WB, IP, IHC-P, IF-F, FC-FP	M
79809T		
Motif		
80353S Carbamidomethyl Cysteine (F3T7J) Rabbit mAb	WB	All
80353T		
95132S β-Hydroxybutyryl Lysine (E6H5Q) Rabbit mAb	W	All
95132T		
24141SF Phospho-Tyrosine (P-Tyr-1000) MultiMab® Rabbit mAb mix (BSA and Azide Free)	WB, IF-IC, FC-FP	All



Learn more at
cellsignal.com/newproducts



Cell Signaling
 TECHNOLOGY®

PRODUCT		APPLICATIONS	REACTIVITY
Neural function (neurite extension, axonal guidance, synaptogenesis, myelination and neuron-glia cell interactions)			
95730S	ALDH1L1 (E712Q) Feline Chimeric mAb	IF-F	H, M, R
95730T			
30803S	ALDH1L1 (E712Q) Horse Chimeric mAb	IF-F	H, M, R
30803T			
49973S	ALDH1L1 (E712Q) Mouse Chimeric mAb	IF-F	H, M, R
49973T			
96510S	BASP1 (F9A5E) Rabbit mAb	WB, IHC-P, IF-IC	H
96510T			
91686S	Phospho-CREB (Ser133) (87G3) Feline Chimeric mAb	IF-F	H, M, R
91686T			
54388S	Phospho-CREB (Ser133) (87G3) Horse Chimeric mAb	IF-F	H, M, R
54388T			
87358S	Phospho-CREB (Ser133) (87G3) Mouse Chimeric mAb	IF-F	H, M, R
87358T			
23180S	Dab1 (F6E1K) Rabbit mAb	WB, IP	H, M, R
23180T			
99059S	Doublecortin (F6K9E) Feline Chimeric mAb	IF-F	H, M, R
99059T			
30298S	Doublecortin (F6K9E) Horse Chimeric mAb	IF-F	H, M, R
30298T			
41762S	Doublecortin (F6K9E) Mouse Chimeric mAb	IF-F	H, M, R
41762T			
35460SF	Neurofilament-H (E9C9Z) Mouse mAb (BSA and Azide Free)	WB, IF-F, IF-IC	H, M, R
18593S	OMG (F8B7X) Rabbit mAb	WB, IHC-P	H
18593T			
67931S	S100B (E7C3A) Rabbit mAb (Alexa Fluor® 488 Conjugate)	IHC-P, IF-IC	H, M, R
76475S	S100B (E7C3A) Rabbit mAb (Alexa Fluor® 647 Conjugate)	IHC-P, IF-F, IF-IC	H, M, R
76103S	TMEM119 (E3E1O) Feline Chimeric mAb	IF-F	M
76103T			
88925S	TMEM119 (E3E1O) Horse Chimeric mAb	IF-F	M
88925T			
26113S	TMEM119 (E3E1O) Mouse Chimeric mAb	IF-F	M
26113T			
92081S	β3-Tubulin (D71G9) XP® Feline Chimeric mAb	IF-F	H, M, R
92081T			
59820S	β3-Tubulin (D71G9) XP® Horse Chimeric mAb	IF-F	H, M, R
59820T			
70668S	β3-Tubulin (D71G9) XP® Mouse Chimeric mAb	IF-F	H, M, R
70668T			
Neurodegeneration: Other			
95895S	HDGF2 (Cryptic Specific) Antibody	WB	H
93535P	Huntingtin (D7F7) XP® Rabbit mAb (Biotinylated)	E-P	H, M, R
60987S	Phospho-TMEM106B (Tyr50) (F3L6S) Rabbit mAb	WB, IP	H, M, R
60987T			
NFκappaB Signaling			
94741S	TRAF5 (F6W4N) Rabbit mAb	WB, IHC-P	H, M
94741T			
Nuclear Structure			
20049SF	PML (E6S9L) Rabbit mAb (BSA and Azide Free)	IHC-P, IF-IC	H
Parkinson's disease			
77511SF	CDK5 (D1F7M) Rabbit mAb (BSA and Azide Free)	WB, IF-F, IF-IC	H, M, R, Mk
Phenotyping and Functional Markers			
50036S	CD31 (PECAM-1) (F8M3S) Rabbit mAb (Alexa Fluor® 488 Conjugate)	IHC-P	H, Mk
70753S	CD31 (PECAM-1) (F8M3S) Rabbit mAb (Alexa Fluor® 555 Conjugate)	IHC-P	H, Mk
20964S	CD31 (PECAM-1) (F8M3S) Rabbit mAb (Alexa Fluor® 647 Conjugate)	IHC-P	H, Mk
15822S	CD34 (QBEND/10) Mouse mAb	IHC-P	H
91388S	MHC Class II (I-Eα) (F8N7J) Rabbit mAb	WB, IP, IHC-P	M
91388T			
83905S	Neutrophil Elastase (E8U3X) Rabbit mAb (Alexa Fluor® 488 Conjugate)	IHC-P	M
71796S	Neutrophil Elastase (E8U3X) Rabbit mAb (Alexa Fluor® 647 Conjugate)	IHC-P	M
PI3K/Akt/mTOR Signaling			
81368SF	SHP-2 (F4V3P) Rabbit mAb (BSA and Azide Free)	WB, IHC-P	H, M, R
73020SF	Phospho-S6 Ribosomal Protein (Ser235/236) (E2R1O) Mouse mAb (BSA and Azide Free)	WB, IF-IC	H, M, R, Mk
Post-synapse			
78512S	Neurexin-3 (F6Q9O) Rabbit mAb	WB, IP	H, M, R
78512T			
Pre-synapse & transmitters			
82985S	Pancreatic polypeptide (F8J4Y) Rabbit mAb	WB, IHC-P, IF-F	H, M
82985T			



Learn more at
cellsignal.com/newproducts



Cell Signaling
 TECHNOLOGY®

PRODUCT	APPLICATIONS	REACTIVITY
Pre-synapse & transmitters (cont.)		
76788S Synaptophysin (D8F6H) XP® Feline Chimeric mAb	IF-F	H, M, R
76788T		
94064S Synaptophysin (D8F6H) XP® Horse Chimeric mAb	IF-F	H, M, R
94064T		
99929S Synaptophysin (D8F6H) XP® Mouse Chimeric mAb	IF-F	H, M, R
99929T		
Proteostasis		
37155P CRBN (D8H3S) Rabbit mAb (Biotinylated)	E-P	H, M, R
46155S ISG15 (F2W5O) Rabbit mAb	WB, IHC-P, FC-FP	H
46155T		
18486S ISG15 (F3L3J) Rabbit mAb	WB	H
18486T		
60857P KEAP1 (D6B12) Rabbit mAb (Biotinylated)	E-P	H, M, R, (Mk, B, Pg)
99483S Prostate Specific Membrane Antigen (D7I8E) XP® Rabbit mAb (Alexa Fluor® 488 Conjugate)	IHC-P	H, M, R
73418S Prostate Specific Membrane Antigen (D7I8E) XP® Rabbit mAb (Alexa Fluor® 647 Conjugate)	IHC-P	H, M, R
67517SF SERPINB3 (F3Z8E) Rabbit mAb (BSA and Azide Free)	WB, IHC-P, SEC	H
56408S SLPI (F5M6K) Rabbit mAb	WB	M
56408T		
32920S UBE3C (F2S2J) Rabbit mAb	WB, IP, IF-IC	H, M, R
32920T		
Receptor signaling: GPCRs, ion channels, RTKs-Neuro, transporters (SLC & ABC)		
56678S GPR75 (F3P5Y) Rabbit mAb	WB	H, M, R
56678T		
63168S Kv1.3/KCNA3 (F2R6N) Rabbit mAb	WB, IP	H, M, R
63168T		
42042S TMEM35 (F3G2T) Rabbit mAb	WB, IP	H, M, R
42042T		
52971S TrkB (F8L9G) Rabbit mAb	WB, IHC-P	H, M, R
52971T		
92447SF TrkC (C44H5) Rabbit mAb (BSA and Azide Free)	WB, IHC-P, IF-IC	H, M, R
REDOX/Energy Metabolism		
97366S Phospho-ATP-Citrate Lyase (Ser455) (F4T8R) Rabbit mAb	WB, IP, IHC-P, FC-FP	H, M, R
97366T		
80448SF DLAT (4A4-B6-C10) Mouse mAb (BSA and Azide Free)	WB, IF-IC	H, Mk
63366SF OPA1 (D7C1A) Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H
RNA Modifications (Editing, Methylation, PseudoU, Other)		
96964T Human Reactive ADAR1 Antibody Sampler Kit	-	-
RNA Processing		
15588SF TIAR (D32D3) XP® Rabbit mAb (BSA and Azide Free)	WB, IF-IC, FC-FP	H, M, R, Mk
RNA Stability & Decay		
83312S SMG5 (F9P6O) Rabbit mAb	WB, IP	H, Mk
83312T		
RTK-Cell Biology		
95087SF Phospho-ALK (Tyr1507) (D6F1V) Rabbit mAb (BSA and Azide Free)	WB, IF-IC, FC-FP	H, (M, R)
94995SF Phospho-FGF Receptor 1 (Tyr653/654) (D4X3D) Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H, (M, R)
TAG		
37497S Myc-Tag (71D10) Feline Chimeric mAb	IF-IC	All
37497T		
81163S Myc-Tag (71D10) Horse Chimeric mAb	IF-IC	All
81163T		
94152S SRLEEELRRRLTE Tag (F8Z6D) Rabbit mAb (Binds to same epitope as NanoTag Biotechnologies Anti-ALFA antibody)	WB, IF-IC	All
94152T		
31431S V5-Tag (D3H8Q) Feline Chimeric mAb	IF-IC	All
31431T		
46509S V5-Tag (D3H8Q) Horse Chimeric mAb	IF-IC	All
46509T		
67168S V5-Tag (D3H8Q) Mouse Chimeric mAb	IF-IC	All
67168T		
23765SF VSV-G Tag (E8S5G) Rabbit mAb (BSA and Azide Free)	WB, IF-IC	All
TGFB/BMP		
94665S Activin Receptor Type 1 (F6M1E) Rabbit mAb	WB, IP	H, M, R, Mk
94665T		
72877SF SMAD2 (86F7) Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H, Mk
Transcription Factors (Hematopoiesis and Leukemogenesis)		
32303S IRF-8 (F7M3K) Rabbit mAb	WB, IP, IHC-P	M, R
32303T		
47603SF POU2F2 (F3B7M) Rabbit mAb (BSA and Azide Free)	WB, IF-F, IF-IC, FC-FP	H, M
37929SF PU.1 (F2D5A) Mouse mAb (BSA and Azide Free)	WB, IF-F, IF-IC, FC-FP	M



Learn more at
cellsignal.com/newproducts



Cell Signaling
 TECHNOLOGY®

PRODUCT	APPLICATIONS	REACTIVITY
Transcription Factors (Hematopoiesis and Leukemogenesis) (cont.)		
94103S ZBTB32 (F501L) Rabbit mAb	WB, FC-FP	M
94103T		
Translational Control		
90731SF Sec61B (D5Q1W) Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H, M, R, Mk
Tumor Microenvironment		
35382S CD47 (D3O7P) Rabbit mAb (Alexa Fluor® 488 Conjugate)	IHC-P	H
24205S CD47 (D3O7P) Rabbit mAb (Alexa Fluor® 647 Conjugate)	IHC-P	H
66825S CXCR4 (FIK2M) Rabbit mAb	WB, IHC-P	M
66825T		
34845S MHC Class II (LGII-612.14) Mouse mAb (Alexa Fluor® 488 Conjugate)	IHC-P, FC-L	H
53814S Mouse IL-10 Recombinant Protein	-	-
69324SF TNF-α (D1G2) Rabbit mAb (BSA and Azide Free)	WB, IF-IC, FC-FP	H
Vesicle and Protein Trafficking		
44700SF Phospho-AP2M1 (Thr156) (D4F3) Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H, (M, R, Mk)
43781S CHMP7 (F4P9J) Rabbit mAb	WB, IP	H
43781T		
38736SF Clathrin Heavy Chain (D3C6) XP® Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H, M, R, Mk, (Pg)
24923SF GM130 (D6B1) XP® Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H, Mk
44573C PathScan® RP NPM1 (C Mutant Specific) Sandwich ELISA Kit	ELISA	H
57398SF Rab10 (D36C4) XP® Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H, M, R, Mk
73006SF Rab11 (D4F5) XP® Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H, M, R, Mk, (Hm, Pg)
84334S Rab32 (F6E5T) Rabbit mAb	WB, IP	M
84334T		
39384SF Transferrin Receptor/CD71 (D7S5Z) Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H
WNT		
88194SF DKK1 (D5V6L) Rabbit mAb (BSA and Azide Free)	WB, IF-IC	H
InTraSeq™ シングルセル解析関連製品		
83302S InTraSeq™ 3' Conjugate Antibody Cocktail 2	SCA-ITS3	H, M
85815S Bcl-xL (54H6) Rabbit mAb (InTraSeq™ 3' Conjugate 3058)	SCA-ITS3	H, M
31103S Cleaved Caspase-3 (Asp175) (D3E9) Rabbit mAb (InTraSeq™ 3' Conjugate 3078)	SCA-ITS3	H, M
17166S c-Fos (E217R) XP® Rabbit mAb (InTraSeq™ 3' Conjugate 3057)	SCA-ITS3	H, M
44436S Phospho-c-Fos (Ser32) (D82C12) XP® Rabbit mAb (InTraSeq™ 3' Conjugate 3056)	SCA-ITS3	H, M
66731S c-Jun (60A8) Rabbit mAb (InTraSeq™ 3' Conjugate 3072)	SCA-ITS3	H, M
81686S Phospho-c-Jun (Ser63) (E617P) XP® Rabbit mAb (InTraSeq™ 3' Conjugate 3070)	SCA-ITS3	H, M
69852S Ezh2 (D2C9) XP® Rabbit mAb (InTraSeq™ 3' Conjugate 3081)	SCA-ITS3	H, M
26379S FGF Receptor 2 (D4L2V) Rabbit mAb (InTraSeq™ 3' Conjugate 3071)	SCA-ITS3	H, M
99585S HIF-1α (D1S7W) XP® Rabbit mAb (InTraSeq™ 3' Conjugate 3066)	SCA-ITS3	H, M
92853S Histone H2A.X (D17A3) XP® Rabbit mAb (InTraSeq™ 3' Conjugate 3082)	SCA-ITS3	H, M
83332S Histone H3 (D1H2) XP® Rabbit mAb (InTraSeq™ 3' Conjugate 3051)	SCA-ITS3	H, M, Hm, (C, Dm, X, Z, B)
32730S Histone H3 (D1H2) XP® Rabbit mAb (InTraSeq™ 3' Conjugate 3052)	SCA-ITS3	H, M, Hm, (C, Dm, X, Z, B)
43286S Acetyl-Histone H3 (Lys9) (C5B11) Rabbit mAb (InTraSeq™ 3' Conjugate 3083)	SCA-ITS3	H, M
65200S Tri-Methyl-Histone H3 (Lys4) (C42D8) Rabbit mAb (InTraSeq™ 3' Conjugate 3085)	SCA-ITS3	H, M
58599S Tri-Methyl-Histone H3 (Lys36) (D5A7) XP® Rabbit mAb (InTraSeq™ 3' Conjugate 3084)	SCA-ITS3	H, M
47881S JunB (C37F9) Rabbit mAb (InTraSeq™ 3' Conjugate 3069)	SCA-ITS3	H, M
99506S Ki-67 (D3B5) Rabbit mAb (InTraSeq™ 3' Conjugate 3076)	SCA-ITS3	H, M
96305S LC3B (E5Q2K) Mouse mAb (InTraSeq™ 3' Conjugate 3062)	SCA-ITS3	H, M
48487S mTOR (7C10) Rabbit mAb (InTraSeq™ 3' Conjugate 3064)	SCA-ITS3	H, M
84197S Phospho-NDRG1 (Thr346) (D98G11) XP® Rabbit mAb (InTraSeq™ 3' Conjugate 3075)	SCA-ITS3	H, M
17838S Notch2 (D76A6) XP® Rabbit mAb (InTraSeq™ 3' Conjugate 3060)	SCA-ITS3	H, M
58425S Phospho-p38 MAPK (Thr180/Tyr182) (D3F9) XP® Rabbit mAb (InTraSeq™ 3' Conjugate 3077)	SCA-ITS3	H, M
83854S Phospho-Rb (Ser807/811) (D20B12) XP® Rabbit mAb (InTraSeq™ 3' Conjugate 3079)	SCA-ITS3	H, M
34360S Phospho-SMAD2 (Ser465/Ser467) (E8F3R) Rabbit mAb (InTraSeq™ 3' Conjugate 3061)	SCA-ITS3	H, M
59131S Phospho-SQSTM1/p62 (Ser349) (E7M1A) Rabbit mAb (InTraSeq™ 3' Conjugate 3074)	SCA-ITS3	H, M
85390S YAP (D8H1X) XP® Rabbit mAb (InTraSeq™ 3' Conjugate 3065)	SCA-ITS3	H, M
62771S Phospho-YAP (Ser61) (E4Q8S) Rabbit mAb (InTraSeq™ 3' Conjugate 3055)	SCA-ITS3	H, M
50265S Whitlow/218 Linker (E3U7Q) Rabbit mAb (InTraSeq™ 3' Conjugate 3053)	SCA-ITS3	All




Learn more at
cellsignal.com/newproducts



Cell Signaling
 TECHNOLOGY®

What's New at CST




The SignalStar® Multiplex IHC Secondary Kit is here!

- Up to 70% faster time-to-results
- Signal amplification of up to 8 markers (even with low expression)
- Automated and flexible panel design
- Now with access to >2,000 IHC-validated antibodies

cst-science.com/Q3-2025-info3

Remarkable **flexibility**.
Now available on demand.



New SignalFlex™ Antibody Conjugates

Signal Flex™ offer you a broader range of fluorochrome-antibody combinations to choose from, using the same primary CST® antibodies you trust. Enjoy greater flexibility when building your fluorescent immunoassay panels for applications like flow cytometry, with:

- More fluorochrome options
- Easy to order
- Convenient size

cst-science.com/Q3-2025-info2

© 2024 Cell Signaling Technology, Inc. CST, Cell Signaling Technology, and InTraSeq are registered trademarks of Cell Signaling Technology, Inc. 10x Genomics, 10x, Feature Barcode, and Chromium are the trademarks or registered trademarks of 10x Genomics, Inc. All other trademarks are the property of their respective owners. Visit cellsignal.com/trademarks for more information. Subject to patents licensed from 10x Genomics, Inc., for use with single-cell (i.e., Chromium) 10x products.



Learn more at
cellsignal.com/newproducts



Cell Signaling
TECHNOLOGY®