


SSI is introducing DuraFrame™ PCR Plates: The Premium Choice for Quality and Value



At the heart of DuraFrame lies **unparalleled durability**, engineered to withstand the rigors of modern laboratories. Our plates feature a rugged **polycarbonate** frame, providing a solid and flat surface that guarantees seamless interaction with automated liquid handling platforms. With a firm grip assured, robots can confidently place DuraFrame PCR Plates with utmost precision, eliminating any concerns of misalignment.

Not stopping at durability alone, we have designed DuraFrame with **polypropylene wells** that ensure **optimal thermal transfer** and **higher reaction efficiency**. This translates to more accurate and reliable results, giving researchers the confidence they need in their experiments. Furthermore, our **white well options** have been meticulously crafted to **enhance the fluorescence signal** during qPCR, enabling clearer and more precise measurements.

Sealing PCR plates is a crucial aspect of laboratory workflows, and DuraFrame excels in this domain as well. The **slightly raised wells** of DuraFrame PCR Plates guarantee an **evaporation-proof seal** with film, foil, or strip caps. No more worries about sample integrity or loss due to evaporation. DuraFrame ensures your valuable samples are protected throughout the entire process.

With an **ANSI/SBS footprint and well spacing**, DuraFrame PCR Plates seamlessly integrate into existing laboratory setups, offering compatibility and ease of use. Moreover, our **vibrant color options** provide a visual distinction between different samples or experiments, allowing for quick and efficient identification and organization.

DuraFrame Introductory Offer:

Please check reverse side for all available plate types

Description	Quantity	List Price CHF	Offer Price CHF
DuraFrame 96-well (all colors frames/clear wells)	10 Plates/Pack	77.00	57.80
DuraFrame 96-well (clear frame/white wells)	10 Plates/Pack	88.50	66.40
DuraFrame 384-well (all colors frames/clear wells)	10 Plates/Pack	106.50	79.90
DuraFrame 384-well (clear frame/white wells)	10 Plates/Pack	110.50	82.90

This promotion is valid until December 31st 2024. No other discounts apply.
All prices in CHF excl. VAT. Please mention offer number **SDF1224** when ordering.

DuraFrame Introductory Offer:

Cat. No.	Description	Quantity	List Price CHF	Offer Price CHF
96-well Full-Skirt 0.1mL PCR Plates				
3742-00	Duraframe 96-well Full-Skirt 0.1mL PCR Plate, clear frame/clear wells	10 plates	77.00	57.80
3742-01	Duraframe 96-well Full-Skirt 0.1mL PCR Plate, blue frame/clear wells	10 plates	77.00	57.80
3742-02	Duraframe 96-well Full-Skirt 0.1mL PCR Plate, green frame/clear wells	10 plates	77.00	57.80
3742-04	Duraframe 96-well Full-Skirt 0.1mL PCR Plate, red frame/clear wells	10 plates	77.00	57.80
3742-05	Duraframe 96-well Full-Skirt 0.1mL PCR Plate, violet frame/clear wells	10 plates	77.00	57.80
3742-06	Duraframe 96-well Full-Skirt 0.1mL PCR Plate, yellow frame/clear wells	10 plates	77.00	57.80
3742-09	Duraframe 96-well Full-Skirt 0.1mL PCR Plate, assorted (B,G,R,V,Y) frames/clear wells	10 plates	77.00	57.80
3742-W0	Duraframe 96-well Full-Skirt 0.1mL PCR Plate, clear frame/white wells	10 plates	88.50	66.40
96-well Semi-Skirt 0.2mL PCR Plates				
3750-00	Duraframe 96-well Semi-Skirt 0.2mL PCR Plate, clear frame/clear wells	10 plates	77.00	57.80
3750-01	Duraframe 96-well Semi-Skirt 0.2mL PCR Plate, blue frame/clear wells	10 plates	77.00	57.80
3750-02	Duraframe 96-well Semi-Skirt 0.2mL PCR Plate, green frame/clear wells	10 plates	77.00	57.80
3750-04	Duraframe 96-well Semi-Skirt 0.2mL PCR Plate, red frame/clear wells	10 plates	77.00	57.80
3750-05	Duraframe 96-well Semi-Skirt 0.2mL PCR Plate, violet frame/clear wells	10 plates	77.00	57.80
3750-06	Duraframe 96-well Semi-Skirt 0.2mL PCR Plate, yellow frame/clear wells	10 plates	77.00	57.80
3750-09	Duraframe 96-well Semi-Skirt 0.2mL PCR Plate, assorted (B,G,R,V,Y) frames/clear wells	10 plates	77.00	57.80
384-well Full-Skirt PCR Plates				
3730-00	Duraframe 384-well Full-Skirt PCR Plate, clear frame/clear wells	10 plates	106.50	79.90
3735-00	Duraframe 384-well Full-Skirt PCR Plate, clear frame/clear wells	10 plates	106.50	79.90
3735-01	Duraframe 384-well Full-Skirt PCR Plate, blue frame/clear wells	10 plates	106.50	79.90
3735-02	Duraframe 384-well Full-Skirt PCR Plate, green frame/clear wells	10 plates	106.50	79.90
3735-04	Duraframe 384-well Full-Skirt PCR Plate, red frame/clear wells	10 plates	106.50	79.90
3735-05	Duraframe 384-well Full-Skirt PCR Plate, violet frame/clear wells	10 plates	106.50	79.90
3735-06	Duraframe 384-well Full-Skirt PCR Plate, yellow frame/clear wells	10 plates	106.50	79.90
3735-09	Duraframe 384-well Full-Skirt PCR Plate, assorted (B,G,R,V,Y) frames/clear wells	10 plates	106.50	79.90
3735-W0	Duraframe 384-well Full-Skirt PCR Plate, clear frame/white wells	10 plates	110.50	82.90

This promotion is valid until December 31st 2024. No other discounts apply.
All prices in CHF excl. VAT. Please mention offer number **SDF1224** when ordering.

Interested in testing? Please contact us for a free sample!

DuraFrame PCR plates are also available **pre-barcoded**:

Cat. No.	Description	Quantity	List Price CHF
3730T00	Duraframe 384-well Full-Skirt PCR Plate, clear frame/clear wells, pre-barcoded	10 plates	115.00
3735T00	Duraframe 384-well Full-Skirt PCR Plate, clear frame/clear wells, pre-barcoded	10 plates	115.00
3735TW0	Duraframe 384-well Full-Skirt PCR Plate, clear frame/white wells, pre-barcoded	10 plates	119.00
3742T00	Duraframe 96-well Full-Skirt 0.1mL PCR Plate, clear frame/clear wells, pre-barcoded	10 plates	85.50
3742TW0	Duraframe 96-well Full-Skirt 0.1mL PCR Plate, clear frame/white wells, pre-barcoded	10 plates	97.00



Paradiesrain 14
CH-4123 Allschwil 1

Postfach
Switzerland

Tel. +41 (0)61 486 80 80
Fax +41 (0)61 486 80 00

info@bioconcept.ch
www.bioconcept.ch